LABNOVATION

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- 1.) Nanshan Chen*, Min Zhou*, Xuan Dong*, Jieming Qu*, Fengyun Gong, Yang Han, Yang Qiu, Jingli Wang, Ying Liu, Yuan Wei, Jia'an Xia, Ting Yu, Xinxin Zhang, Li Zhang Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study. LANCET. January 29, 2020.
- 2.) World Health Organization (Coronavirus disease 2019) [https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-\(covid-2019\)-and-the-virus-that-causes-it\(Zugriff am 27.03.2020\)](https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-(covid-2019)-and-the-virus-that-causes-it(Zugriff am 27.03.2020))
- 3.) World Health Organization (Coronavirus disease 2019) https://www.who.int/health-topics/coronavirus#tab=tab_1 (Zugriff am 27.03.2020)

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INSTRUCTION FOR USE

INTRODUCTION

In December 2019, the novel respiratory disease (COVID-19) caused by the coronavirus (SARS-CoV-2) was reported in Wuhan, China.^{1,2} According to WHO, most of the people infected with SARS-CoV-2 have mild to moderate respiratory diseases, fever, cough and recover without treatment. However, people with weak immune systems, such as the elderly or people with previous illnesses (e.g., cardiovascular disease, diabetes, chronic respiratory diseases, cancer, etc.) are more likely to develop a serious illness that can lead to the death of the infected person.³

This rapid test kit is intended for the qualitative detection of SARS-CoV-2 viral nucleocapsid antigens from human anterior nasal of secretion from individuals suspected of COVID-19. Positive result of the antigen test can be used for early isolation of patients with suspected infection, but it cannot be used as diagnosis basis of SARS-CoV-2 infection. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment. Further nucleic acid detection should be carried out for suspected population whose antigen test result is positive or negative.

This kit is an immunochromatography assay which detects SARS-CoV-2 nucleocapsid antigen in the samples with the help of the double antibody sandwich method. If there is virus antigen presence in the sample, it binds with the corresponding colloidal gold antibody. This complex "migrates" across the membrane and binds to the monoclonal antibody at the Test line (T). This creates a visible red line, which indicates a positive result. However, if the sample does not contain any antigen, then the complex cannot be formed and thus no reddish line forms in the Test line (T). Regardless of whether the sample contains antigen or not, a reddish line forms in the Control line (C).

KIT COMPONENTS

- > 1 Test cassette
- > 1 Sample tube with prefilled sample extraction buffer
- > 1 Swab
- > 1 Tube stand (in the outer box)
- > 1 Instruction for use

Additionally required materials:

1 Timer

STORAGE & SHELF LIFE

The product should be stored at 2-30°C, and its shelf life is 24 months. The test cassette and sample extraction buffer should be used immediately after opening.

EVALUATION OF TEST RESULTS

To read the test results simply determine whether a line is present or absent at the Control (C) position. It does not matter how strong or weak a Control line (C) is.

PERFORMANCE CHARACTERISTICS

1. Sensitivity and Specificity

Analytical results for all samples from clinical studies:

| SARS-CoV-2 Antigen Rapid Test | RT-PCR | | TOTAL |
|-------------------------------|----------|----------|-------|
| | Positive | Negative | |
| Positive | 306 | 0 | 306 |
| Negative | 8 | 650 | 658 |
| TOTAL | 314 | 650 | 964 |

Sensitivity: > 95%*

Specificity: >99 %*

Total accuracy: >98% *

*95 % Confidence Interval.

2. Limit of detection

| | |
|-------------------|---------------------------|
| LOD concentration | 30 TCID ₅₀ /mL |
|-------------------|---------------------------|

3. Cross-reactivity

With human coronavirus 229E, human coronavirus OC43, human

coronavirus NL63, adenovirus, human metapneumovirus, parainfluenza virus 1, parainfluenza virus 2, parainfluenza virus 3, parainfluenza virus 4, influenza A, type B Influenza, enterovirus, respiratory syncytial virus, rhinovirus, SARS coronavirus, MERS coronavirus, Haemophilus influenzae, Streptococcus pneumoniae, Streptococcus pyogenes, Candida albicans, Bordetella pertussis, Mycoplasma pneumonia, pneumonia Chlamydia, Legionella pneumophila, etc. have no cross reaction.

4. Interfering

Common interfering substances in the sample, such as blood, mucin, and pus, have no effect on the test results.

WARNINGS AND IMPORTANT INFORMATION

- This kit is a qualitative detection, which cannot determine the exact content of antigen.
- The test is intended for use outside the body only.
- Not to be taken internally. Avoid sample buffer contact with skin and eyes.
- Protect from sunlight, do not freeze. Store in a dry place between 2°C and 30°C. Do not use after the expiration date printed on the package.
- Keep out of the reach of children. Any child under age 16 shouldn't perform the test without parental guidance, or professional aid.
- Not following the exact instructions can affect the outcome of the test. The final diagnosis must be confirmed by a physician.
- Do not use the test if the packaging is damaged. Do not use broken test components.
- All test components are only intended to be used for this test. Do not reuse the test or test components.
- The test should be carried out immediately or within one hour after opening the foil pouch (15-30°C, humidity <60%).
- Samples be processed as soon as possible after sample collection. If the test cannot be performed immediately, the sample should be stored in a sealed state, stored at 2~8°C for 8 hours, and stored below -20°C for 1 month. Long-term storage is not recommended.
- Poor vision, color blindness or poor lighting may affect your ability to interpret the test correctly.
- DISPOSAL The test kit can be disposed of with normal household waste in accordance with applicable local regulations.
- A negative result does not rule out the infection of a SARS-CoV-2 infection. Therefore, the test should not be used as the only reference for the clinical diagnosis. The result must be confirmed by the PCR.
- After use, rinse hands or, in case of contact with the buffer solution, the affected body parts thoroughly with water.
- If symptoms persist: Seek medical advice.

LITERATURE

1.) Nanshan Chen*, Min Zhou*, Xuan Dong*, Jieming Qu*, Fengyong Gong, Yang Han, Yang Qiu, Jingli Wang, Ying Liu, Yuan Wei, Jia'an Xia, Ting Yu, Xinlin Zhang, Li Zhang Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study. LANCET. January 29, 2020.

2.) World Health Organization (Coronavirus disease 2019) [https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-\(covid-2019\)-and-the-virus-that-causes-it](https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-(covid-2019)-and-the-virus-that-causes-it) (Zugriff am 27.03.2020)

3.) World Health Organization (Coronavirus disease 2019) https://www.who.int/health-topics/coronavirus#tab=tab_1 (Zugriff am 27.03.2020)

STEP-BY-STEP INSTRUCTION

Choose a location to do this test where it can sit UNDISTURBED for 20 minutes.

Let test cassette and test components stand at a room temperature (10°C to 30°C) before performing the test.

① Wash and dry hands before you begin to perform the test.

② Please check the expiration date printed on the BOX Do not use it beyond the expiration date.

Lay all the supplied materials on a clean, dry and flat surface.

③ Take sample tube with sample extraction buffer and remove the sealed film of the sample tube. Then place the tube in the tube stand (The tube stand is the outer box, see below.)

④ Open the sealed pouch and remove the test cassette. Lay it face up on a clean, dry and flat surface.

⑤ Remove the swab from the container, being careful NOT to touch the soft end, which is the absorbent tip.


⑥ Gently, insert the entire absorbent tip of the swab (around 1.5 cm) into your nostril. Slowly, rotate the swab in a circular against the inside walls of your nostril 5 times or more. Be sure to collect any nasal drainage that maybe present on the swab. Gently remove the swab. Use the same swab to repeat steps in the other nostril and slowly, take out the swab.

CAUTION: If the swab stick breaks during specimen collection, repeat specimen collection with a new swab.

⑦ Insert the swab into the sample tube with prefilled extraction buffer. Mix well and squeeze the swab 10-15 times by compressing the walls of the tube against the swab.


Roll the swab head against the inner wall of the tubes as you remove it. Try to release as much liquid as possible. Dispose of the used swab in accordance with applicable local regulations.

⑧ Close the cap of the sample tube.

Add 3 full drops of the mixed solution vertically into the sample well (S) of the test cassette.

Read the result 15-20 minutes after adding the sample. Result got after 20 minutes is invalid.


| | |
|--|--|
| | Keep out of the reach of children. Any child under age 16 shouldn't perform the test without parental guidance, or professional aid. |
| | Improper handling of specimen can lead to biological infections. Avoid any direct contact with reagents or wastes. |

INTERPRETATION OF RESULTS

Positive

Two colored bands appear on the membrane. One band appears in the control region (C) and another band appears in the test region (T). The test result means that SARS-CoV-2 antigen is detectable in your sample. The detection of these antigens indicate with a high probability of infection with

the novel coronavirus.

In case of a positive test result:

- There is currently a suspicion of a COVID-19 infection.
- Immediately contact a doctor/family physician or the local public health department
- follow local guidelines for self-isolation
- have a PCR confirmatory test performed
- *Note: The thickness of the line is insignificant; Any reddish color in the Test line (T) should be considered a positive result. The positive test result must be confirmed by PCR.

Negative

Only one colored band appears, in the control region (C). No apparent colored band appears in the test region (T). The test result indicates that there is no or too little SARS-CoV-2 Antigen in the sample and at the current time there is probably no infection with the novel coronavirus.

Negative

If the test result is negative

- Continue to follow all applicable rules regarding contact with others and protective measures.
- Even if the test is negative, an infection may still be present.
- In case of suspicion, repeat the test after 1 - 2 days, as the coronavirus cannot be detected accurately in all phases of an infection.
- *Note: False negative results can be from incorrect sampling, incorrect execution of the test, or insufficient virus in the sample.

Invalid

If there is no Control line (C) or only a Test line (T) in the result window, the test did not run correctly and the results are not valid.

In case of an invalid test result

- Possibly caused by incorrect test performance
- Repeat the test

- If test results remain invalid, contact a doctor or a COVID-19 testing centre.

*Note: It is important that you carefully follow the instructions for the test. You should test again with a new sample and a new test.

DISPOSE USED TEST

DISPOSAL The test kit can be disposed of with normal household waste in accordance with applicable local regulations.

*Note: After completing all steps, wash hands or use hand sanitizer.

MANUFACTURER

LABNOVATION TECHNOLOGIES, INC.

Address: 101 and 5th Floor, Building 1, No. 68, 18th Road, Guangming Hi-Tech Park, Tangjia Community, Fenghuang Street, Guangming District, Shenzhen 518107, China
Tel: 0086-755-86368398 Web: www.labnovation.com
E-mail: export@labnovation.com

Hotline: +49 251 3226669
Borkstrasse 10, 48163 Münster, Germany

INSTRUCTIONS OF SYMBOL

| | |
|--|--|
| | |
| | |
| | |
| | |
| | |
| | |

Ref Ver.1.0
2021-12-30